


Rapid Rhino devices feature high volume low pressure tamponade to ensure gentle and even compression to control epistaxis at the source of the bleed.

Rapid Rhino has a unique design that will control epistaxis the first time and is more comfortable for patients upon insertion and removal than other packing options. Now Rapid Rhino can control posterior epistaxis with RR 900.



Introducing Rapid Rhino® 900 for posterior epistaxis


Unilateral Models




Rapid Rhino 450
 • 4.5 cm inflatable epistaxis device



Rapid Rhino 750
 • 7.5 cm inflatable epistaxis device



Rapid Rhino 551
 • 5.5 cm inflatable epistaxis device that contains an airway




Rapid Rhino 900 NEW!
 • 9 cm tamponade device for posterior epistaxis
 • Contains 2 balloons to control intranasal hemorrhage originating from the sphenopalatine artery



Rapid Rhino 550
 • 5.5 cm inflatable epistaxis device



Rapid Rhino 530
 • Rapid-Pac™ is 5.5 cm
 • Offers a non-inflatable option for epistaxis




Rapid Rhino 751
 • 7.5 cm inflatable epistaxis device that contains an airway

Bilateral Models



Rapid Rhino 552
 • 5.5 cm bilateral epistaxis device
Rapid Rhino 555
 • 5.5 cm bilateral epistaxis device that contains an airway



Rapid Rhino 752
 • 7.5 cm bilateral epistaxis device
Rapid Rhino 755
 • 7.5 cm bilateral epistaxis device that contains an airway



ArthroCare ENT
 680 Vaqueros Avenue
 Sunnyvale, CA 94085-3523
 main phone 408-736-0224
 order entry **800-797-6520**
 order entry fax 888-994-2782
www.arthrocareENT.com

ArthroCare Europe AB
 Skeppsbron 2
 111 30 Stockholm
 Sweden
 Tel +46 8 546 172 00
 Fax +46 8 546 172 39
info@arthrocare.se

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

© 2009 ArthroCare Corporation. ArthroCare, Rapid-Pac and Rapid Rhino are registered trademarks of ArthroCare Corporation. Rhino Rocket is a registered trademark of Shippert Medical Technologies Corp.

Rapid Rhino is the epistaxis solution

Rapid Rhino is a nasal tamponade device that is specifically designed to address the major challenges in the management of epistaxis. Now introducing Rapid Rhino 900. This 9 cm pack is the first Rapid Rhino device with 2 balloons designed to control posterior bleeds originating from the Sphenopalatine artery. Finally there is a device that controls posterior epistaxis on the spot, while still maintaining the level of patient comfort Rapid Rhino is known for.

Rapid Rhino 900 product usage directions

Soak for minimum 30 seconds




1

Material handling

- To use Rapid Rhino 900, remove the device from the envelope package.
- Remove blue plastic tube encasing.
- Thoroughly saturate the product in *sterile water* by submerging the product for at least a full 30 seconds.

Rapid Rhino device placement




2

Placement technique

- Insert the Rapid Rhino 900 into the patient's nostril parallel to the septal floor, or following along the superior aspect of the hard palate, until the blue indicator ring is inside the opening of the nostril.

Cuff inflation



3


Posterior cuff inflation

- Locate the inflation line with the *green striped* swallow guard butterfly. This indicates the posterior cuff which should be inflated first.
- Using a 20 cc syringe, slowly inflate the posterior balloon with *air only* inside the patient's nose. The balloon will conform to the anatomy of the posterior nasal cavity.
- Use the pilot cuff as you would on an endotracheal tube to monitor intranasal pressure as you inflate.
- Stop inflation when pilot cuff is rounded and firm.

Technique Pearl

When there is sufficient pressure in the posterior chamber, remove the syringe from the inflation valve and gently pull the device anteriorly, as necessary, to ensure there is sufficient pressure on the posterior bleed site.

Gentle and even pressure to the site of the bleed.




4

Anterior cuff inflation

- Next, locate the second inflation line, without a swallow guard butterfly. This is to inflate the anterior balloon. Again, use the pilot cuff to maintain constant tactile feedback of intranasal pressure taking care not to overinflate.

Patient discharge




5

Patient discharge

- When there is sufficient pressure in both balloons allow the patient to sit for 15-20 minutes prior to discharge. Swelling in the nasal anatomy will reduce and the balloons may need to be inflated more to avoid movement of the device.
- After this second assessment of the pilot cuff tape both inflation catheters to the patient's cheek for discharge.

Deflation and removal



6

Removal

- The patient may come back 24-72 hours later to remove the device or they may see an ENT specialist for follow-up.